

## NATtrol™ Respiratory Panel (RSP) Controls

**PRODUCT DESCRIPTION:**

**NATtrol™ Respiratory Panel (RSP) Controls\* (qualitative)** is formulated with purified, intact bacterial cells and viral particles. The microorganisms have been chemically modified to render them non-infectious and refrigerator stable. NATRSPC-BIO contains 6 x 0.3 mL vials of RSP Control 1 and 6 x 0.3 mL of RSP Control 2. See Table 1 for control contents. The controls are supplied in a proprietary matrix.

\*Pat.: <http://www.zeptometrix.com/patent-information/>

**INTENDED USE:**

- NATtrol™ RSP Controls are designed to evaluate the performance of nucleic acid tests for determination of the presence of bacterial and viral nucleic acids (from organisms listed in Table 1). NATtrol™ Respiratory Panel (RSP) Controls enable laboratories to monitor test variation, lot-to-lot test kit performance, operator variation, and can provide assistance in identifying random or systemic error.

**WARNINGS AND PRECAUTIONS:**

- NATtrol™ inactivation was carried out on microorganism stocks used to formulate the control members. The inactivation was verified in a standard microbiological growth protocol.
- This product contains inactivated microorganisms and materials of human and animal origin. Safe practices suggest that the controls be considered potentially infectious and to use Universal Precautions when handling.
- Refer to CDC guidelines and local regulations for handling and disposal.
- The matrix used in the manufacture of this product is treated with 0.09% sodium azide. It was manufactured from Human Serum Albumin that have been tested and found to be non-reactive at the donor level for HIV-1/HIV-2 Antibody, HBsAg and HCV Antibody by FDA licensed donor screening test methods. All materials are also tested for HIV-1 and HCV by FDA approved Nucleic Acid Test (NAT) methods.
- Heat inactivated Fetal Bovine Serum used in the manufacture of this product meet applicable USDA requirements for abattoir sourced animals, traceability and country of origin. The materials were collected at USDA licensed establishments or legally imported from countries recognized by the USDA as negligible or controlled for risk for Bovine Spongiform Encephalopathy (BSE) and other exotic disease agents. Donor animals were inspected ante and post mortem at the abattoir as required by the USDA.
- Do not use past the expiration date on the label.
- To avoid cross-contamination, use separate pipette tips for all materials.

**RECOMMENDED STORAGE:**

- NATtrol™ Respiratory Panel (RSP) Controls should be stored at 2-8°C.

**INSTRUCTIONS FOR USE:**

- Mix vial vigorously for at least 5 secs.
- Process according to manufacturer's instructions for sample to result assays.
- Extract nucleic acid prior to use in downstream assays that are not sample to result.

**LIMITATION:**

- FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**
- Quality control materials should be used in accordance with local, state, federal, and accreditation requirements.
- This product is not intended to replace the manufacturer's controls provided with the assay.

**EXPECTED RESULTS:**

- Each laboratory must evaluate the product and establish their own acceptance criteria.
- These controls have been tested with the BIOFIRE® SPOTFIRE® Respiratory Panel assay on the BIOFIRE® SPOTFIRE® System and provide all expected results for the targets listed in Table 1.
- Table 1 is for informational purposes only.

PINATRSPC-BIO  
Revision: 00  
Effective Date: 03/24/2023

	Catalog Number		Temperature Limitation
	Batch Code		Expiration Date
	For Research Use Only		Biological Risk
	Manufacturer		

PCA# N/A  
Page 1 of 2

# ZeptoMetrix®

## NATtrol™ Respiratory Panel (RSP) Controls

Catalog Number:  
NATRSPC-BIO

**TABLE 1: Organisms/Viruses in Respiratory Panel (RSP) Controls**

Organism/Virus	Strain	RSP Control 1	RSP Control 2
Adenovirus Type 1	N/A	Positive	Negative
Adenovirus Type 3	N/A	Positive	Negative
Adenovirus Type 31	N/A	Positive	Negative
<i>C. pneumoniae</i>	Z500	Positive	Negative
Influenza A H1	A/New Caledonia/20/99	Positive	Negative
Influenza A 2009 H1N1pdm	A/NY/02/2009 <sup>1</sup>	Positive	Negative
Influenza A H3	A/Brisbane/10/07	Positive	Negative
<i>M. pneumoniae</i>	M129	Positive	Negative
Metapneumovirus 8	Peru6-2003 <sup>2</sup>	Positive	Negative
Rhinovirus 1A	N/A	Positive	Negative
SARS-CoV-2	USA-WA1/2020 <sup>3</sup>	Positive	Negative
<i>B. parapertussis</i>	A747	Negative	Positive
<i>B. pertussis</i>	A639	Negative	Positive
Coronavirus 229E	N/A	Negative	Positive
Coronavirus HKU-1	Recombinant <sup>4</sup>	Negative	Positive
Coronavirus NL63	N/A	Negative	Positive
Coronavirus OC43	N/A	Negative	Positive
Influenza B	B/Florida/02/06	Negative	Positive
Parainfluenza Type 1	N/A	Negative	Positive
Parainfluenza Type 2	N/A	Negative	Positive
Parainfluenza Type 3	N/A	Negative	Positive
Parainfluenza Type 4	N/A	Negative	Positive
RSV A	N/A	Negative	Positive

<sup>1</sup> Please note that although similar in nomenclature, this is a 2009 H1N1 pandemic Influenza strain and does NOT correlate with the seasonal 2009 Influenza strains found in the FluDb.org database. For reference, the NCBI Taxon IDs for the seasonal Influenza strains listed in the FluDb.org database are: A/New York/01/2009 (H1N1) - 666252; B/New York/01/2009 - 664512; A/New York/02/2009 (H1N1) - 666298; and A/New York/03/2009 (H3N2) - 659637.

<sup>2</sup> This product is sold by ZeptoMetrix under license from Vironovative B. V. under patent applications, including U.S. Patent Applications 10/371,099 and 10/371,12 and any patents that issue from applications related to PCT/NL02/00040 and PCT/US03/05271.

<sup>3</sup> This reagent was deposited by the Centers for Disease Control and Prevention and obtained through BEI Resources, NIAID, NIH: SARS-Related Coronavirus 2, Isolate USA-WA1/2020, NR-52281.

<sup>4</sup> This analyte only contains a short sequence of the viral genome therefore each laboratory must evaluate performance in their assay.

PINATRSPC-BIO  
Revision: 00  
Effective Date: 03/24/2023

	Catalog Number		Temperature Limitation
	Batch Code		Expiration Date
	For Research Use Only		Biological Risk
	Manufacturer		

PCA# N/A  
Page 2 of 2