

Parainfluenza Virus
Type: 4A
Purified Viral Lysate (1 mg)

Catalog Number: 0810060

### PRODUCT DESCRIPTION:

Parainfluenza Virus (Type: 4A) is an enveloped, negative-sense single-stranded RNA virus with a diameter of approximately 150 to 250 nm.

The virus is purified using sucrose density gradient ultracentrifugation, disrupted in the presence of 0.5% Triton X-100 non-ionic detergent/0.6M KCl, and heat inactivated. Each frozen aliquot contains 1 mg of viral lysate per the BCA total protein assay.

Viral inactivation is verified by the absence of viral growth in tissue culture-based infectivity assays.

### **INTENDED USE:**

Viral lysates are intended for research, product development, or further manufacturing use. The performance characteristics of this product should be tested by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

### **BIOSAFETY:**

Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

# PRECAUTIONS:

- Use Universal Precautions, this product is potentially biohazardous.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary).
- To avoid cross-contamination, use separate pipette tips for all reagents.

## RECOMMENDED STORAGE:

Viral lysates should be stored at -65°C or below.

PI0810060 Revision: 08 Effective Date: 09/13/2022

ı	REF	Catalog Number	1	Temperature Limitation
	LOT	Batch Code	M	Expiration Date
	RUO	For Research Use Only	<b>®</b>	Biological Risk
ſ	-	Manufacturer		•