LeptoMetrix® Human Herpes Virus Type 6B (HHV-6B) Strain: Z29 Infectious Culture Fluid (1 mL)

PRODUCT DESCRIPTION:

Human Herpes Virus Type 6B (HHV-6B) (Strain: Z29) is an enveloped icosahedral virus 150 to 200 nm in diameter that contains a nucleocapsid and double-stranded linear DNA.

Each frozen aliquot contains 1 mL of titered viral culture fluid.

HHV-6 does not always demonstrate distinct CPE in culture. Therefore, HHV-6 culture fluid is assayed by either Tissue Culture Infective Dose (TCID50), or real-time PCR (copies/mL). Please indicate if you prefer TCID50 or PCR quantitation when placing an order.

INTENDED USE:

Viral culture fluids are sold as consumable testing materials; propagation or commercialization is prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).

BIOSAFETY:

HHV-6B is a Biosafety Level 2 microorganism and must be used within a Biological Safety Level 2 facility or cabinet. Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous.**
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer will be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

PI0810072CF Revision: 06 Effective Date: 09/15/2022

| REF | Catalog Number | X | Temperature Limitation |
|-----|-----------------------|---|------------------------|
| LOT | Batch Code | M | Expiration Date |
| RUO | For Research Use Only | 9 | Biological Risk |
| | Manufacturer | | |

PCA# 22-308 Page 1 of 1

www.ZeptoMetrix.com ZeptoMetrix LLC • 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487 This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

Catalog Number: 0810072CF