Adenovirus Type 41 (Species F) Strain: Tak Infectious Culture Fluid (1 mL) Catalog Number: 0810085CF

PRODUCT DESCRIPTION:

Adenovirus Type 41 (Species F) (Strain: Tak) is a non-enveloped icosahedral virus that contains a nucleocapsid and double-stranded linear DNA.

Each frozen aliquot contains 1 mL of titered viral culture fluid.

INTENDED USE:

Viral culture fluids are sold as consumable testing materials, propagation or commercialization are prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).

BIOSAFETY:

Adenovirus is a Biosafety Level 2 microorganism and must be used within a Biological Safety Level 2 faciliy or cabinet. Please consult your institution's regulations regarding the use of this organism. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

PRECAUTIONS:

- Use Universal Precautions, this product is potentially biohazardous.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary).
 Titer will be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

PI0810085CF Revision: 06 Effective Date: 09/09/2022

	REF	Catalog Number	4	Temperature Limitation
	LOT	Batch Code	M	Expiration Date
	RUO	For Research Use Only	€	Biological Risk
	444	Manufacturer		